



InSinkErator **LEDefense**

FLG-SS-LED

Submittal Sheet

LEDefense Features

- Kills 99%¹ of odor-causing bacteria
- UV-Free LED Light
- Always on 24/7 Protection
- Helps prevent buildup on bottom of baffle

Warranty


3-Year Limited Warranty.²

Certifications

Meets UL430.

Specifications

Includes	Sink flange, baffle, power adapter and mounting components.
Retrofit	LEDefense is compatible with all InSinkErator® continuous feed disposals.
Install Fitment	Fits all standard sinks 3/4" thick or less. Cannot be used with sink thickness greater than 3/4". Replaces current sink flange.
Power	LEDefense is powered by a USB-C power adapter with a 6ft cord. LEDs are powered 24/7. The power draw is <4.7W. Requires a continuously powered outlet.
LED Light	LEDefense is UV-free
Product Weight	0.4lbs
Shipping Weight	1.7lbs



Job Specifications



InSinkErator® LEDefense

Directions for LEDefense Use:
Must be continuously powered on for antimicrobial action. If powered off, antimicrobial action stops.

Keep Out of Reach of Children | EPA Est No. 105182-WI-1

Sold by InSinkErator, LLC.
1250 International Dr, Mount Pleasant, WI USA, 53177

1-800-558-5700
www.insinkerator.com

¹ Grind Chamber/Flange Below Baffle: Results may vary depending on the amount of light that is reaching the surfaces in the space where the product is installed and the length of time of exposure. Use of Vyv antimicrobial light is not intended to replace manual cleaning and disinfection practices. Sanitization will only function when unit is plugged in and LED lights are on. Testing on a non-enveloped virus (MS2 bacteriophage) showed a 99.985% reduction in controlled laboratory testing in 6 hours on hard surfaces. Testing on SARS-CoV-2 (enveloped virus) showed a 98.45% reduction in controlled laboratory testing in 4 hours on hard surfaces. Testing on Staphylococcus aureus, Salmonella enterica and Escherichia coli showed 99%+ reduction in controlled laboratory testing in 24 hours on hard surfaces. Growth of mold (Aspergillus) inhibited (in laboratory-controlled studies with prolonged use).

² Visit www.insinkerator.com for warranty details.